



BASISPLATTE / BASE PLATE

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PICTROGRAM LEGEND

en

Pictogrammes listed below you'll find on the product packaging, if applicable.

Catalogue number

Batch code

Use by

Content (in pieces)

Consult instructions for use

Keep away from sunlight

Keep dry

Do not reuse

Medical devices

CE marking

Manufacturer

BASE PLATE

1. FOREWORD

These instructions for use are valid for all base plates.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

2. SAFETY INFORMATION

Base plates are a single-patient product and is intended for single use only.

Observe the respective manufacturer's instructions when using tracheostoma valves.

Only aids with a 22mm adapter may be used, otherwise there is a risk of aspiration.

When used, check compatibility and a secure fit with utmost care.

3. INTENDED PURPOSE

Base plates are an attachment system for daily use of HME filter cassettes or tracheostoma valves with a 22 mm adapter.

4. INDICATION

The product is indicated for tracheotomized and laryngectomized patients.

5. CONTRAINDICATION

Use of the tape should be avoided while undergoing radiotherapy due to the risk of skin irritations.

After undergoing radiotherapy, only use base plates after consulting your treating physician.

Start using the tape at the earliest 2 weeks after completing radiotherapy.

Do not use if the patient is allergic to the material.

6. COMPLICATIONS AND ADVERSE REACTIONS

The circular connector opening can under certain circumstances reduce the usable diameter of the tracheostoma, resulting in a reduction of air supply.

The particular properties of the adhesive used for the medical tape can in individual cases cause skin irritations.

Should you experience skin irritation (e.g. redness, itchiness), discontinue the product immediately and consult a physician.

This may be due to special skin sensitivity or allergy-related incompatibility.

7. PRODUCT DESCRIPTION

Base plates are an attachment system for daily use of HME filter cassettes or tracheostoma valves. The self-adhesive tapes are selected according to the anatomical situation of the tracheostoma and the type of skin in order to achieve a long wearing period and optimal attachment.

The filter cassette is fitted in place in the **22mm** combined adapter and can be changed separately from the base plate.

8. APPLICATION

Check whether the packaging and the product are in a faultless condition.

If a product is damaged, it must be discarded immediately.

Check the use-by or expiry date.

Do not use the product after this date.

1. Clean the skin around the tracheostoma in front of a mirror using warm water and pH-neutral soap or **OPTIFAHL®** cleaning cloths (**REF 33200**).

Dry the skin.

Before cleaning, remove residues of adhesive from the skin if necessary, using **OPTICLEAR®** (**REF 33500**).

2. To protect the skin from irritation, apply **OPTIGARD®** (**REF 33600**) to the affected area of the skin.
3. The base plates should be warmed up for **3** minutes before use by holding them between the palms of the hands in order to optimise the adhesive properties and to facilitate adaptation to the affected area of skin (see picture **1**).
4. After removing the protective foil from the reverse side, (see picture **2**) bring the lower edge of the opening of the base plates into position level with the lower edge of your tracheostoma (see picture **3**) and then fix the tape in place on the skin.
5. Smooth this area with your finger to ensure that the base plate is positioned correctly.
6. Now remove the other section of the protective film (if such exists) and smooth the area with your finger (see picture **4**).
7. To ensure optimal fixation and to avoid air enclosures, carefully stroke the entire tape surface with your fingers for **5** minutes after attaching the tape to the skin.
Optimal skin attachment is achieved after 20 minutes, only start inserting the filter cassette when this time has elapsed.
This prevents unwanted air enclosures between the base plate and the skin.
8. To remove a filter cassette, fix the base plates in **CONVEX** position with **2** fingers (see picture **5**) and carefully take out the filter cassette in order to achieve a long wearing period of the tape.
9. To remove the base plates, carefully grip an adhesive tab.
Carefully peel this off from one side, if necessary with **OPTICLEAR®** (**REF 33500**) (see picture **6**).

CAUTION!

Do not use creams that contain fat! Otherwise, secure attachment of the tape is no longer ensured.

9. HYGIENE INSTRUCTIONS

Cleaning, disinfection, re-sterilisation and re-use can impair the safety and function of the product and are therefore not permissible!

10. STORAGE

Store in a dry environment away from sunlight and/or heat.

11. SERVICE LIFE

The maximum period of use / wearing period is **24 hours**.

When skin irritations occur, the base plate must be changed immediately.

Deviations from the expected wearing period occur as a result of individual wearing conditions and skin properties.

12. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

13. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section **11**, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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