



Laryngotec Pro Kombi Lingo tube maat 12/36 - 14722-1236

<https://mdsmedical.nl/laryngotec-pro-kombi-lingo-tube-maat-1236-14722-1236>

Specificaties

Productnummer	14722-1236
Categorie	Buttons & Tubes
Model	Laryngotec Pro Kombi Lingo
Kleur	transparant
Materiaal	silicone
Type	gefenestreerd
Verpakkingseenheid	per stuk
Merk	FAHL
Laryngectomie Tube/Button Type	Laryngotec Pro Kombi Lingo
Maat	12-36

LARYNGOTEC® PRO KOMBI LINGO: Laryngectomie tube maat 12/36